



TORSIONAL VIBRATION APPARATUS

HVT13



Year 3
study

Features

- Very visual demonstration of Torsional vibration
- Fully guarded motor
- Free and fixed end conditions
- Digital oscilloscope and Software supplied as standard
- Manual Twist of specimen using set of weights and hanger
- Mechanical Damping supplied
- Oscillation speed digitally controlled and displayed

Description

This bench top mounted unit is used to study torsion and torsional vibration. The apparatus is constructed around a profiled aluminium base with levelling feet, onto which is mounted 4 vertical chuck pillars. Each pillar contains a central shaft running in precision bearings with a chuck at one end which grips a torsion specimen. Each chuck also contains a large disc of varying mass and inertia. With these chucks it is possible to assemble a torsion specimen with up to 3 masses. Vibrations are transmitted into the torsion specimen by means of an exciter, which is electronically speed controlled from the main control unit. To change the end conditions of the apparatus a fifth chuck houses a chuck which rigidly clamps the end of the torsion specimen to achieve a fixed end. Also supplied with the apparatus is a manual torsion arrangement, which allows a known angular twist to be applied to the specimen. A cord is wrapped around one of the large discs and a load is applied via a hanger and weights set. The angle of twist for incremental loading is recorded and the modulus of rigidity can be calculated. Oscillations sensors are mounted integrally with each mass pillar and provide signals of the amplitude of vibration. The control unit conditions these signals and makes them available to the supplied digital oscilloscope for vibration analysis.

Related Laws/Applications

- Torsion
- Torsional stiffness
- Resonance
- Multi-mass torsional systems

Learning capabilities

- Determination of the torsional stiffness of a torsion bar
- Determination of the mass moment of inertia Decay behaviour of torsional vibration
- Forced torsional vibration and resonance
- Torsional vibration systems with three weights allowing single, two and three mass torsional vibration to be investigated.
- Demonstration of the effect of frictional damping

Technical Specification

- Stainless steel torsion specimen: 1300mm long x Ø6mm
- Crank offset: 0.15mm (approximately)
- Up to 20Hz excitation frequency
- 20mm thick x Ø150mm inertia disc (2 off), Steel
- 15mm thick x Ø228mm inertia disc (1 off), plastic
- 6 x 5N Weights

What's in the Box?

- 1 x HVT13
- 1 x Control console
- 1 x Power lead
- 5 x Spare fuse
- 4 x BNC cable
- 1 x Spare timing belt
- 1 x Tape measure
- 1 x Stop watch
- 1 x Load hanger
- 1 x Hex Wrench set
- Instruction manual
- Packing list
- Test sheet
- Digital Oscilloscope

You might also like

- HSM2
- HSM31
- HVT12b

Weights & Dimensions

- Weight: 55 kg
- Length: 1400mm
- Width: 410mm
- Height: 400mm

Essential Services

- Sturdy bench top
- 110/120V, 60Hz or 220/240V, 50Hz, single phase, live neutral and earth

Ordering information

To order this product, please call PA Hilton quoting the following code: HVT13

All brand and/or product names are trademarks of their respective owners. Specifications and external appearance are subject to change without notice. The colour of the actual product may vary from the colour shown in the brochure.

Copyright © 2018 P.A. Hilton Limited. All rights reserved. This technical leaflet, its contents and/or layout may not be modified and/or adapted, copied in part or in whole and/or incorporated into other works without the prior written permission of P. A. Hilton Limited. Hi-Tech Education is a registered trade mark of P. A. Hilton Limited.

COUNTRY OF ORIGIN - UK WARRANTY PERIOD - 5 YEARS